



DEVICE DETAILS

NAME OF DEVICE	WRIST DIGITAL BLOOD PRESSURE MONITOR
ESTABLISHMENT NAME	REVO HEALTHCARE SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GB7946820-40403
BRAND NAME	A&D
MEDICAL DEVICE CATEGORY	MD 1302 - Monitoring devices of vital physiological parameters
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	The A&D Wrist Digital Blood Pressure Monitors use an oscillometric method, and detects pulse rate and determines blood pressure by pressurizing a cuff around the wrist. When the start button is pressed, the device detects pulse rate during inflation, and decides pressurization setting value, and automatically inflates up to the pressurization value. If the pre-set pressurization is not enough to detect patient's BP value, the device automatically repeats re-inflation continuously. The device also has constant deflation, automatic quick deflation and other safety functions to make it easier to operate.
DEVICE INTENDED PURPOSE	The A&D Wrist Digital Blood Pressure Monitor is intended to measure the systolic, diastolic and pulse rate of an individual by using noninvasive technique in which an inflatable cuff is wrapped around the wrist.
VALIDITY DATE OF REGISTRATION	25/02/2020 - 24/02/2025

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	UB-525 Wrist Digital Blood Pressure Monitor	UB-525
2	UB-533 Wrist Digital Blood Pressure Monitor	UB-533
3	UB-533PGMR Wrist Digital Blood Pressure Monitor	UB-533PGMR