



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	A&D UPPER ARM BLOOD PRESSURE MONITOR
<b>ESTABLISHMENT NAME</b>	REVO HEALTHCARE SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	AUTHORIZED REPRESENTATIVE
<b>REGISTRATION NO</b>	GB2631319-34850
<b>BRAND NAME</b>	A&D
<b>MEDICAL DEVICE CATEGORY</b>	MD 1302 - Monitoring devices of vital physiological parameters
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	The A&D Blood Pressure Monitor is a device designed to measure blood pressure and pulse rate using non-invasive technique in which an inflatable cuff is wrapped around the user's upper arm. The technique to define the blood pressure (systolic and diastolic pressures) is similar to the auscultatory method but utilize an electronic capacitive pressure sensor rather than the mercury manometer and stethoscope. This sensor converts the tiny alteration in the cuff pressure to electrical signals; these signals are analyzed to define the systolic and diastolic pressures and calculating pulse rate which is a well known technique called the oscillometric method.
<b>DEVICE INTENDED PURPOSE</b>	The A&D Blood Pressure Monitor is intended to measure the systolic, diastolic and pulse rate of an individual by using noninvasive technique in which an inflatable cuff is wrapped around the upper arm.
<b>VALIDITY DATE OF REGISTRATION</b>	19/09/2019 - 18/09/2024

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	UA-611 BLOOD PRESSURE MONITOR	UA-611
2	UA-651 BLOOD PRESSURE MONITOR	UA-651
3	UA-651SL BLOOD PRESSURE MONITOR	UA-651SL
4	UA-767JP BLOOD PRESSURE MONITOR	UA-767JP
5	UA-767F BLOOD PRESSURE MONITOR	UA-767F
6	UA-767F-W BLOOD PRESSURE MONITOR	UA-767FW
7	ADULT CUFF	CUF-F-A
8	LARGE ADULT CUFF	CUF-F-LA
9	SEMI LARGE ADULT CUFF	CUF-F-SL

<b>NO</b>	<b>NAME OF DEVICE</b>	<b>IDENTIFIER</b>
10	WIDE RANGE CUFF	CUF-I
11	UA-611 BLOOD PRESSURE MONITOR	UA-611
12	ADULT CUFF	CUF-F-A
13	LARGE ADULT CUFF	CUF-F-LA
14	UA-611 UPPER ARM BLOOD PRESSURE MONITOR	UA-611
15	ADULT CUFF	CUF-F-A
16	LARGE ADULT CUFF	CUF-F-LA
17	UA-651 UPPER ARM BLOOD PRESSURE MONITOR	UA-651
18	ADULT CUFF	CUF-F-A
19	LARGE ADULT CUFF	CUF-F-LA
20	UA-651SL UPPER ARM BLOOD PRESSURE MONITOR	UA-651SL
21	SEMI LARGE ADULT CUFF	CUF-F-SL
22	LARGE ADULT CUFF	CUF-F-LA
23	UA-767JP UPPER ARM BLOOD PRESSURE MONITOR	UA-767JP
24	ADULT CUFF	CUF-F-A
25	LARGE ADULT CUFF	CUF-F-LA
26	WIDE RANGE CUFF	CUF-I
27	UA-767F UPPER ARM BLOOD PRESSURE MONITOR	UA-767F
28	ADULT CUFF	CUF-F-A
29	LARGE ADULT CUFF	CUF-F-LA
30	UA-767F-W UPPER ARM BLOOD PRESSURE MONITOR	UA-767F-W
31	ADULT CUFF	CUF-F-A
32	LARGE ADULT CUFF	CUF-F-LA
33	WIDE RANGE CUFF	CUF-I