



DEVICE DETAILS

NAME OF DEVICE	PARTICULATE RESPIRATOR KN95
ESTABLISHMENT NAME	PENMED MARKETING SDN. BHD.
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA5294221-75140
BRAND NAME	ZHUSHI
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	mainly consists of the body of the mask, assemble the mask belt and nose clip. the outer layer of the mask is PP non-woven cloth, the middle layer is white hot air cotton, double layer filter melt spray cloth, hot press to make.
DEVICE INTENDED PURPOSE	Used to protect the wearer from airborne particles and splashes of potentially contaminated liquids which intended for surgical, medical and general use.
VALIDITY DATE OF REGISTRATION	08/10/2021 - 07/10/2026

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	PARTICULATE RESPIRATOR KN95	ZS9501