



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	KN95 MEDICAL FACE MASK
<b>ESTABLISHMENT NAME</b>	ONEWORLD CORPORATE (M) SDN. BHD.
<b>ROLE OF ESTABLISHMENT</b>	AUTHORIZED REPRESENTATIVE
<b>REGISTRATION NO</b>	GA7339621-75836
<b>BRAND NAME</b>	SURGIPLUS
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	SINGLE
<b>DEVICE DESCRIPTION</b>	KN95 Medical face mask is to be use by healthcare professional during procedures to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids, and particulate material. - BFE >95% FFP2 5 -layers surgical & protective mask - Barrier layer: Non-woven fabric blocking large particles - Filtration layer:Melt-blown layer adsorb small particles - Filtration layer: Electrostatic capture of PM0.3 particles
<b>DEVICE INTENDED PURPOSE</b>	This KN95 Medical Face Mask effectively filters to serve as barrier to potential air pollutants that can be harmful to patient /surgeon and healthcare personnel. The face mask's filtration efficiency of 5 layers is intended to prevent air pollutants like allergens, virus, bacteria, dust, droplets and other harmful non-oil based particles in the air. The Face mask is Suitable to be used in Surgical Room , Medical Clinic and Hospitals
<b>VALIDITY DATE OF REGISTRATION</b>	13/10/2021 - 12/10/2026

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	KN95 MEDICAL FACE MASK	KN95 Mask (5ply)