



DEVICE DETAILS

NAME OF DEVICE	KN 95 FACE MASK
ESTABLISHMENT NAME	GSH PRECISION TECHNOLOGY SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA8786121-76250
BRAND NAME	KN95
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	COMES WITH EXTERNAL LAYER POLYPROPYLENE MATERIAL WATER REPELLENT LAYER AND INTERNAL LAYER NON WOVEN AND EP FABRIC MATERIAL
DEVICE INTENDED PURPOSE	GSH-KN95 Medical Mask is intended to be worn by health professionals during surgery and certain health care procedures to catch microorganisms shed in liquid droplets and aerosols from the wearer's mouth and nose. The medical personnel are dressed with them in operation so that the blood, body fluid and exudation of patients, which has potential infection, can be obstructed effectively to prevent the bacteria infection. All parts of the products and processing materials are used to meet the requirements of medical materials, convenient use, safety, etc.
VALIDITY DATE OF REGISTRATION	18/10/2021 - 17/10/2026

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	KN95/FFP2 Medical Face Mask	GSH KN95/FFP2 FM 001
2	KN95/FFP2 Medical Face Mask	GSH KN95/FFP2 FM 002
3	KN95/FFP2 Medical Face Mask	GSH KN95/FFP2 FM 003
4	KN95/FFP2 Medical Face Mask	GSH KN95/FFP2 FM 004
5	KN95/FFP2 Medical Face Mask	GSH KN95/FFP2 FM 005
6	KN95/FFP2 Medical Face Mask	GSH KN95/FFP2 FM 006
7	KN95/FFP2 Medical Face Mask	GSH KN95/FFP2 FM 007
8	KN95/FFP2 Medical Face Mask	GSH KN95/FFP2 FM 008
9	KN95/FFP2 Medical Face Mask	GSH KN95/FFP2 FM 009

NO	NAME OF DEVICE	IDENTIFIER
10	KN95/FFP2 Medical Face Mask	GSH KN95/FFP2 FM 010
11	KN95/FFP2 Medical Face Mask	GSH KN95/FFP2 FM 011



MEDICAL DEVICE AUTHORITY