

## DEVICE DETAILS

NAME OF DEVICE	Clefmed Premium KN95 Surgical Face Mask
ESTABLISHMENT NAME	CLEFTECH TRADING SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA4912822-94424
BRAND NAME	CLEFMED+ (HYNAUT)
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	Foldable, 4-ply with ear loops and equipped with a minimum filter efficiency of BFE, PFE & VFE 99%. Composed of high-quality materials such as non-woven, melt-blown fabric.
DEVICE INTENDED PURPOSE	As a facial protective mask which helps to shield against particulate matter inclusive of germs and bacteria, and any exposure towards pollutant in air including air-borne droplets through mask's physical barriers.
VALIDITY DATE OF REGISTRATION	13/06/2022 - 12/06/2027

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Clefmed Premium KN95 Surgical Face Mask	C005