



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	NEUTROVIS KN95 FACE RESPIRATOR
<b>ESTABLISHMENT NAME</b>	NEUTROVIS SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	MANUFACTURER
<b>REGISTRATION NO</b>	GA7413323-116670
<b>BRAND NAME</b>	NEUTROVIS
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	Neutrovis KN95 Face Respirator made with 2 layers of inner melt-blown fabric with $\geq 98\%$ BFE & $\geq 98\%$ PFE to able effectively filter viable and nonviable particles such as bacteria, pollen and dust particles form our surrounding
<b>DEVICE INTENDED PURPOSE</b>	Neutrovis KN95 Face Respirator is made to limit unfiltered air entry to provide better protection against airborne particles. The tight fit around the bridge of your nose and sides of your cheeks directs air to pass through the filter so that you're safe to breathe
<b>VALIDITY DATE OF REGISTRATION</b>	24/01/2023 - 23/01/2028

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Neutrovis KN95 Face Respirator	NTV03
2	Neutrovis KN95 Face Respirator	NTV134
3	Neutrovis KN95 Face Respirator	NTV135
4	Neutrovis KN95 Face Respirator	NTV136
5	Neutrovis KN95 Face Respirator	NTV137
6	Neutrovis KN95 Face Respirator	NTV249