



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	NEUTROVIS KN95 FACE RESPIRATOR WITH NOSE PAD
<b>ESTABLISHMENT NAME</b>	NEUTROVIS SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	MANUFACTURER
<b>REGISTRATION NO</b>	GA1675423-116726
<b>BRAND NAME</b>	NEUTROVIS
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	Neutrovis KN95 Face Respirator is made to limit unfiltered air entry to provide better protection against airborne particles. its contoured shape reduces friction between the respirator and mouth for better communication and breathability. The added fitted nose cushion creates a tight seal and does not cause your glasses to fog up
<b>DEVICE INTENDED PURPOSE</b>	Neutrovis KN95 Face Respirator made with 2 layers of inner melt-blown fabric with $\geq 98\%$ BFE & $\geq 98\%$ PFE to able effectively filter viable and nonviable particles such as bacteria, pollen and dust particles form our surrounding
<b>VALIDITY DATE OF REGISTRATION</b>	24/01/2023 - 23/01/2028

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Neutrovis KN95 Face Respirator with Nose Pad	NTV138
2	Neutrovis KN95 Face Respirator with Nose Pad	NTV139
3	Neutrovis KN95 Face Respirator with Nose Pad	NTV140
4	Neutrovis KN95 Face Respirator with Nose Pad	NTV182
5	Neutrovis KN95 Face Respirator with Nose Pad	NTV183
6	Neutrovis KN95 Face Respirator with Nose Pad	NTV175
7	Neutrovis KN95 Face Respirator with Nose Pad	NTV176
8	Neutrovis KN95 Face Respirator with Nose Pad	NTV177
9	Neutrovis KN95 Face Respirator with Nose Pad	NTV178
10	Neutrovis KN95 Face Respirator with Nose Pad	NTV179

<b>NO</b>	<b>NAME OF DEVICE</b>	<b>IDENTIFIER</b>
11	Neutrovis KN95 Face Respirator with Nose Pad	NTV180
12	Neutrovis KN95 Face Respirator with Nose Pad	NTV181



**MEDICAL DEVICE AUTHORITY**