



DEVICE DETAILS

NAME OF DEVICE	KN95 RESPIRATOR
ESTABLISHMENT NAME	WATSON'S PERSONAL CARE STORES SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA7347723-141366
BRAND NAME	WATSONS
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	A single use, non sterile KN95 respirator made of spunbonded polypropylene and meltdown polypropylene. Each KN95 respirator measuring 155mm x 110mm and they are supplied in a pack of 10 pieces, individually pack.
DEVICE INTENDED PURPOSE	The disposable folding mask is intended to be worn to protect personnel from the transfer of particulate material. This mask is suitable for use in protection against the non-toxic solid and liquid aerosols. The length of time this mask can be used depends on contaminants present. This mask may be used for performing patient procedures, or when patients are in isolation (clean environment, sterile cares, processing department, EC and ICU for bedside procedures, etc.) to reduce the risk of spread of infections. As it is comparable to FFP2 grade mask, it can be used for protection against SARS-Covid 19 as well.
VALIDITY DATE OF REGISTRATION	20/07/2023 - 19/07/2028

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	KN95 RESPIRATOR	22370